

SURGICAL MILK


Instrument lubricant and rust inhibitor

Concentrate

GET THE MOST FROM YOUR INSTRUMENTS



2 L 4 L



Effective

Bacteriostatic.

Prevents spotting, staining and instrument rust. Keeps instruments shiny and prevents dulling.

Effective in steam, dry or chemical vapour sterilizers.

Routine use reduces cost of instrument repair and replacement.

Clinically tested to be steam penetrable and steam sterilizable. Will not interfere with the microbial killing efficacy of autoclaving.

Not sticky, oily or greasy.

Convenient

1 part solution to 5 parts water can be used for up to 7 days.

Can be used on a variety of metals: steel, carbide, tungsten, stainless steel and aluminum.

Use During Autoclaving on:

- Dental/Surgical Instruments
- Hinges
- Box Locks
- General Equipment

4 L bottle of concentrate replaces 20 L of ready-to-use solution



Untreated



Treated with

SURGICAL MILK

RECOMMENDED

by

BMT

Chirurgische Instrumente

www.BMTsurgical.com

Global provider of
Surgical Instruments.



For more information please contact your District Sales Manager or call Germiphene 1.800.265.9931 | www.germiphene.com